

## Quick User Guide




**1** Table Of Contents

- 1. Specifications.....1
- 2. Physical description.....2
- 3. IMPORTANT SAFEGUARDS.....3
- 4. Quick User Guide.....4
- 5. Operation Guide.....6
- 6. Important safety precautions before use...8
- 7. What's in the box?.....11
- 8. Frequently Asked Questions.....14
- 9. Maintenance.....15
- 10. Troubleshooting.....16

### Specifications

Size: 5" x 3.8" x 3.1" (128mm x 96mm x 80mm)  
Weight: 0.73 pound (330 grams with batteries)  
Max speed: 1.7"/sec (4.5cm/sec)  
Wheels: TPU  
Fingerettes: Nylon  
Warranty: one year  
Safety regulation: CE  
Batteries: 3 x AA (1.5V)

**2** Physical description



- 1 Power & Mode button
- 2 Top finger motor
- 3 Detachable wheels
- 4 Battery compartment
- 5 Sensors (Reserved)
- 6 Top finger base
- 7 Top Finger nails

**5** Quick User Guide

Select the desired program using the applicable power/mode button. Lay down horizontally on your stomach as flat as possible.



Place WheeMe on your back. WheeMe will start massaging and steering itself over your back.

WheeMe will start the massage sequence and automatically turns itself off after 15 minutes.

WheeMe can be turned off by pressing and holding the Power/mode button for a short moment. Indications: WheeMe's red LED flashes rapidly for one second and then LED is turned off. All motors are turned off.

**3** Physical description



### IMPORTANT SAFEGUARDS:

This product is a non-professional appliance designed to provide gentle massage to the skin. This appliance is designed for household use only. When using an electrical appliance, basic precautions should always be Followed. READ ALL INSTRUCTIONS BEFORE USING

**6** Operation Guide

1. Pressing the program button a fourth time returns WheeMe to the first program (Vibration massage).
2. Every press of the power/mode button starts a New 15 minute session.
3. WheeMe enters into pause mode if it is positioned upside down. WheeMe exits the pause mode when positioned face up again.



4. WheeMe automatically shuts itself off if batteries are too low.

**4** Quick User Guide

Insert the 3 x AA batteries into WheeMe's bottom battery compartment.



Pressing the power button: powers WheeMe up. Once WheeMe is on, one of three programs can be selected by pressing the power button:

1. **Vibrating massage** without using the top fingernails, the default mode. Indications: as WheeMe vibrates, a red LED rapidly flashes.
2. **Relaxation** using WheeMe's top fingernails (Fingernails (7) need to be mounted to finger base (6) first). Best feeling is reached while not wearing a shirt or blouse. Indications: as WheeMe's top fingernails spin, a red LED slowly flashes.
3. **Dancing**. No other massage elements are used. Best feeling is reached while not wearing a shirt or blouse. Indications: as WheeMe's wheels spin, the light of the red LED remains steady.

**7** Operation Guide

5. The WheeMe has a power-saving feature for the batteries when WheeMe is turned off.
6. Trying to power up WheeMe when the batteries are too low causes the red LED to flash for a short moment. The WheeMe then turns itself off.
7. While WheeMe is working, it is possible to interactively force WheeMe to climb up to your neck area for grabbing it and holding it with your hand. In order to do so, raise your back slowly and wait for a few seconds. Grab WheeMe with your hand when it reaches your neck area.



8. Please note that in order to feel WheeMe's top finger and wheels as they move along your body, you should not wear a shirt or blouse.

## 8 Important safety precautions before use !

1. It doesn't happen frequently but occasionally long hair may become entangled in the wheels. Before using WheeMe, tie your hair up or move it to the side of your head .



2. If not used according to instructions, WheeMe might fall off your body. For more information, read the FAQs and visit [www.wheeme.com](http://www.wheeme.com).

3. Don't allow WheeMe to drop on to the floor; it's apt to break. When using WheeMe it is advisable to lie on a carpet, blankets or other soft materials.

## 11 Important safety precautions before use !

17. When travelling or going from one place to another, the device should be stored without batteries in order to prevent it from operating.

18. Keep WheeMe away from direct sunlight

19. The product could scratch or wound the skin of the user and users should be careful and take appropriate protective measures, particularly when the device is being used near the face or sensitive body organs.

20. The company is not responsible for any bodily harm as a result of improper or inappropriate use.

### What's in the box?

1. WheeMe v2.0.
2. Detachable spinning top fingernails.
3. Quick User Guide.

## 14 Frequently Asked Questions

### 3. The top finger tickles me too much.

If you are ticklish, we recommend purchasing the WheeMe's "Feather Fingernails" which are made of softer material. Check availability of the "Feather Fingernails" and other accessories Under our online store: [www.wheeme.com](http://www.wheeme.com)

### 4. I don't enjoy the top finger.

For full pleasure, the top finger should be used when you're not wearing a shirt or blouse.

## 9 Important safety precautions before use !

4. The WheeMe fingernails are fragile. Take extra care when using and storing it.

5. The WheeMe can only detect the curves on a human's back. It's not designed to operate on hard surfaces or to detect the edge of a table.

6. Close supervision is necessary when this appliance is used on or near children or disabled persons.

7. WheeMe contains small parts. Keep away from children under age 3.

8. WheeMe is not intended for sick individuals or babies.

9. Keep away from water at all time.

10. WheeMe is not intended for medical use. Nor is it a medical instrument or professional appliance. It is designed to provide gentle and calming massage to the skin. Do not use this product as a substitute for medical attention.

## 12 Frequently Asked Questions

### 1. Why does WheeMe keep falling off my body?

For best performance you should be lying flat as possible.

### 2. Why is WheeMe's travel distance so limited? Why doesn't it travel all over my back?

For best performance you should be lying flat as possible.

To extend WheeMe's range over your whole back, select the desired program using the power/mode button together with holding WheeMe at a vertical position.



## 15 Maintenance

Maintaining your WheeMe is easy. Simply release the wheels and finger from the main WheeMe body. Clean the WheeMe gently with a moistened towel, taking care to remove any hair attached to the axis between the wheels and the device. Attach the parts back to the main WheeMe device and you're ready to go.



## 10 Important safety precautions before use !

11. It is recommended to use WheeMe 15 minutes at a time and to let the device cool off between sessions.

12. WheeMe should NEVER be used by an individual suffering from any physical ailment that would limit the user's capacity to operate the controls.


13. The use of WheeMe is the responsibility of the user alone. It should be used appropriately and according to instructions.

14. Please keep in mind that WheeMe could possibly become entangled in clothes or delicate materials the user may be wearing.

15. The device should not be washed in water or wetted in any way. It should only be cleaned according to the guidelines on the packaging.

## 13 Frequently Asked Questions

To return back to normal WheeMe's range over your back, select the desired program using the power/mode button together with holding WheeMe at a horizontal position. Alternatively, power off and on WheeMe again. Normal WheeMe's range is restored. Please note that in this adventure mode, WheeMe travels at a higher angle of tilt which might cause it to fall off your body from time to time.



## 16 Troubleshooting

### 1. WheeMe does not turn on or keeps powering down after a short period of use.

Replace the three AA batteries. For best performance and longer hours of fun, don't use cheap batteries.

For heavy users, we recommend AA 1.2V rechargeable batteries.

### 2. WheeMe keeps falling off my back when using according to the instructions.

Replace the three AA batteries. For best performance and longer hours of fun, don't use cheap batteries.

### 3. My hair gets tangled in the wheels and/or the fingernails.

Shut down WheeMe and release the wheels by pulling out the required wheel.

For more trouble shooting and technical support, please visit our website: [www.WheeMe.com](http://www.WheeMe.com)